

PACO34710

Product Information

Size:

50ug

Reactivity:

Agrobacterium fabrum

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in DNA transformation; mediates the nuclear uptake of single-stranded DNA copies of the transferred DNA (T-DNA) element. Binds single-stranded but not double-stranded DNA regardless of nucleotide sequence composition.

Gene ID:

virE2

Uniprot

P08062

Synonyms:

Single-strand DNA-binding protein (63.5 kDa virulence protein), virE2

Immunogen:

Recombinant Agrobacterium fabrum Single-strand DNA-binding protein (1-556AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A