

PACO34734

---

## Product Information

**Size:**

50ug

**Reactivity:**

Staphylococcus aureus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Staphylococcal enterotoxins cause the intoxication staphylococcal food poisoning syndrome. The illness is characterized by high fever, hypotension, diarrhea, shock, and in some cases death.

**Gene ID:**

entC3

**Uniprot**

P0A0L5

**Synonyms:**

Enterotoxin type C-3 (SEC3), entC3

**Immunogen:**

Recombinant Staphylococcus aureus Enterotoxin type C-3 protein (28-266AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A