

# entG Antibody



PACO34738

---

## Product Information

**Size:**

50ug

**Reactivity:**

Staphylococcus aureus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Staphylococcal enterotoxins cause the intoxication staphylococcal food poisoning syndrome. The illness is characterized by high fever, hypotension, diarrhea, shock, and in some cases death.

**Gene ID:**

entG

**Uniprot**

P0A0L7

**Synonyms:**

Enterotoxin type G (SEG), entG, seg

**Immunogen:**

Recombinant Staphylococcus aureus Enterotoxin type G protein (26-258AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A