

PACO34746

Product Information

Size:

50ug

Reactivity:

Bacillus thuringiensis subsp. Tenebrionis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Promotes colloid-osmotic lysis by binding to the midgut epithelial cells of Coleoptera.

Gene ID:

cry3Aa

Uniprot

P0A379

Synonyms:

Pesticidal crystal protein Cry3Aa (73 kDa crystal protein) (Crystalline entomocidal protoxin) (Insecticidal delta-endotoxin CryIIIa(a)), cry3Aa, bt13 cry3A cryC cryIIIa cryIIIa(a)

Immunogen:

Recombinant Bacillus thuringiensis subsp. Tenebrionis Pesticidal crystal protein cry3Aa protein (58-644AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A