

hup1 Antibody



PACO34750

Product Information

Size:

50ug

Reactivity:

Streptomyces lividans

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Histone-like DNA-binding protein which is capable of wrapping DNA to stabilize it, and thus to prevent its denaturation under extreme environmental conditions.

Gene ID:

hup1

Uniprot

P0A3H6

Synonyms:

DNA-binding protein HU 1 (HSI), hup1, hup

Immunogen:

Recombinant Streptomyces lividans DNA-binding protein HU 1 protein (1-93AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A