

mpt63 Antibody



PACO34754

Product Information

Size:

50ug

Reactivity:

Mycobacterium tuberculosis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:1000-1:5000

Protein Background:

Gene ID:

mpt63

Uniprot

P9WIP1

Synonyms:

Immunogenic protein MPT63 (16 kDa immunoprotective extracellular protein) (Antigen MPT63), mpt63


Immunogen:

Recombinant Mycobacterium tuberculosis Immunogenic protein MPT63 protein (30-159AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Western blot. All lanes: mpt63 antibody at 2 μ g/ml. Lane 1: Recombinant Mycobacterium tuberculosis Immunogenic protein MPT63 protein 1 μ g. Lane 2: Recombinant Mycobacterium tuberculosis Immunogenic protein MPT63 protein 0.1 μ g. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 17 kDa. Observed band size: 17 kDa.