

pth Antibody, Biotin conjugated



PACO34769

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

The natural substrate for this enzyme may be peptidyl-tRNAs which drop off the ribosome during protein synthesis. Involved in lambda inhibition of host protein synthesis. PTH activity may, directly or indirectly, be the target for lambda bar RNA leading to rap cell death.

Gene ID:

pth

Uniprot

P0A7D1

Synonyms:

Peptidyl-tRNA hydrolase (PTH) (EC 3.1.1.29), pth

Immunogen:

Recombinant Escherichia coli Peptidyl-tRNA hydrolase protein (1-194AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A