

PACO34774

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play a role in DNA repair. It seems to be involved in an RecBC-independent recombinational process of DNA repair. It may act with RecF and RecO.

Gene ID:

recR

Uniprot

P0A7H6

Synonyms:

Recombination protein RecR, recR

Immunogen:

Recombinant Escherichia coli Recombination protein RecR protein (1-201AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A