

PACO34798

---

## Product Information

**Size:**

50ug

**Reactivity:**

Escherichia coli

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Involved in the assembly of lipopolysaccharide (LPS). Required for the translocation of LPS from the inner membrane to the outer membrane. May form a bridge between the inner membrane and the outer membrane, via interactions with LptC and LptD, thereby facilitating LPS transfer across the periplasm.

**Gene ID:**

IptA

**Uniprot**

P0ADV1

**Synonyms:**

Lipopolysaccharide export system protein LptA, IptA, yhbN

**Immunogen:**

Recombinant Escherichia coli Lipopolysaccharide export system protein LptA protein (28-185AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A