

PACO34802

Product Information

Size:

50ug

Reactivity:

Enterobacter aerogenes

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Ribitol dehydrogenase has been purified to homogeneity from several strains of *Klebsiella aerogenes*.

Gene ID:

rbtD

Uniprot

P00335

Synonyms:

Ribitol 2-dehydrogenase (RDH) (EC 1.1.1.56), rbtD

Immunogen:

Recombinant *Enterobacter aerogenes* Ribitol 2-dehydrogenase protein (1-249AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A