

PACO34822

Product Information

Size:

50ug

Reactivity:

Staphylococcus aureus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

This protein is probably a specific topoisomerase involved in initiating replication. This protein is specifically required and may be rate-limiting for replication of the plasmid in vivo.

Gene ID:

repD

Uniprot

P03065

Synonyms:

Replication initiation protein, repD

Immunogen:

Recombinant Staphylococcus aureus Replication initiation protein (1-311AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A