

NP Antibody



PACO34830

Product Information

Size:

50ug

Reactivity:

Influenza B virus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

NP

Uniprot

P04665

Synonyms:

Nucleoprotein (Nucleocapsid protein) (Protein N), NP

Immunogen:

Recombinant Influenza B virus Nucleoprotein protein (1-560AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A