

Mup11 Antibody



PACO34838

Product Information

Size:

50ug

Reactivity:

Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Binds pheromones that are released from drying urine of males. These pheromones affect the sexual behavior of females.

Gene ID:

Mup11

Uniprot

P04938

Synonyms:

Major urinary protein 11, Mup11, Mup9

Immunogen:

Recombinant Mouse Major urinary proteins 11 and 8 protein (1-151AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A