

PACO34842

Product Information

Size:

50ug

Reactivity:

Human herpesvirus 1

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Chemokine-binding protein that inhibits neutrophils' chemotaxis.

Gene ID:

gG

Uniprot

P06484

Synonyms:

Envelope glycoprotein G (gG) (gG-1), gG

Immunogen:

Recombinant Human herpesvirus 1 Envelope glycoprotein G protein (25-189AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A