

SNRPB2 Antibody



PACO34870

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Involved in pre-mRNA splicing as component of the spliceosome. Associated with sn-RNP U2, where it contributes to the binding of stem loop IV of U2 snRNA.

Gene ID:

SNRPB2

Uniprot

P08579

Synonyms:

U2 small nuclear ribonucleoprotein B" (U2 snRNP B"), SNRPB2


Immunogen:

Recombinant Human U2 small nuclear ribonucleoprotein B" protein (1-225AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human prostate cancer using PACO34870 at dilution of 1:100.