

SNRPB2 Antibody, FITC conjugated



PACO34872

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Involved in pre-mRNA splicing as component of the spliceosome. Associated with sn-RNP U2, where it contributes to the binding of stem loop IV of U2 snRNA.

Gene ID:

SNRPB2

Uniprot

P08579

Synonyms:

U2 small nuclear ribonucleoprotein B" (U2 snRNP B"), SNRPB2

Immunogen:

Recombinant Human U2 small nuclear ribonucleoprotein B" protein (1-225AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A