

PACO34890

---

## Product Information

**Size:**

50ug

**Reactivity:**

Bacillus thuringiensis subsp. kurstaki

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Promotes colloid osmotic lysis by binding to the midgut epithelial cells of many lepidopteran larvae.

**Gene ID:**

cry1Ab

**Uniprot**

P0A370

**Synonyms:**

Pesticidal crystal protein Cry1Ab (130 kDa crystal protein) (Crystalline entomocidal protoxin) (Insecticidal delta-endotoxin CryIA(b)), cry1Ab, bt2 cry-1-2 cry1A(b) cryIA(b) cryIC1

**Immunogen:**

Recombinant Bacillus thuringiensis subsp. kurstaki Pesticidal crystal protein Cry1Ab protein (1022-1155AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A