

PACO34898

Product Information

Size:

50ug

Reactivity:

Mycobacterium tuberculosis

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

The ESAT-6 proteins are secreted by M. tuberculosis and play important roles in pathogenesis.

Gene ID:

esxB

Uniprot

P9WVK4

Synonyms:

ESAT-6-like protein EsxB (10 kDa culture filtrate antigen CFP-10) (Secreted antigenic protein MTSA-10), esxB, cfp10 lhp mtsA10

Immunogen:

Recombinant Mycobacterium tuberculosis ESAT-6-like protein EsxB protein (2-100AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A