

PACO34906

---

## Product Information

**Size:**

50ug

**Reactivity:**

Escherichia coli

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

ATP, dATP, and GTP are equally effective as phosphate donors. CMP and dCMP are the best phosphate acceptors.

**Gene ID:**

cmk

**Uniprot**

P0A6I0

**Synonyms:**

Cytidylate kinase (CK) (EC 2.7.4.25) (Cytidine monophosphate kinase) (CMP kinase) (Protein MssA) (p25), cmk, mssA ycaF ycaG

**Immunogen:**

Recombinant Escherichia coli Cytidylate kinase protein (1-227AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A