

PACO34910

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

cytosol, glutamate-cysteine ligase activity, metal ion binding, cellular response to arsenic-containing substance, cellular response to mercury ion, glutathione biosynthetic process, hyperosmotic response.

Gene ID:

gshA

Uniprot

P0A6W9

Synonyms:

Glutamate--cysteine ligase (EC 6.3.2.2) (Gamma-ECS) (GCS) (Gamma-glutamylcysteine synthetase), gshA, gsh-I

Immunogen:

Recombinant Escherichia coli Glutamate--cysteine ligase protein (1-518AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A