

PACO34914

Product Information

Size:

50ug

Reactivity:

Escherichia coli

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Transcriptional repressor of the nikABCDE operon. Is active in the presence of excessive concentrations of intracellular nickel.

Gene ID:

nikR

Uniprot

P0A6Z6

Synonyms:

Nickel-responsive regulator, nikR, yhhG

Immunogen:

Recombinant Escherichia coli Nickel-responsive regulator protein (1-133AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A