

IFNA1 Antibody, HRP conjugated



PACO34943

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Produced by macrophages, IFN-alpha have antiviral activities. Interferon stimulates the production of two enzymes: a protein kinase and an oligoadenylate synthetase.

Gene ID:

IFNA1

Uniprot

P01562

Synonyms:

Interferon alpha-1/13 (IFN-alpha-1/13) (Interferon alpha-D) (LeIF D), IFNA1; IFNA13

Immunogen:

Recombinant Human Interferon alpha-1/13 protein (24-189AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A