

PACO34970

---

## Product Information

**Size:**

50ug

**Reactivity:**

Escherichia coli

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

This methylase recognizes the double-stranded sequence GATATC, causes specific methylation on A-2 on both strands, and protects the DNA from cleavage by the EcoRV endonuclease.

**Gene ID:**

ecoRVM

**Uniprot**

P04393

**Synonyms:**

Modification methylase EcoRV (M. EcoRV) (EC 2.1.1.72) (Adenine-specific methyltransferase EcoRV), ecoRVM

**Immunogen:**

Recombinant Escherichia coli Modification methylase EcoRV protein (1-298AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A