

PACO34978

Product Information

Size:

50ug

Reactivity:

Enterobacteria phage T4

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Plays a role in viral DNA synthesis by promoting enzymatic activities of UvsX recombinase, by promoting UvsX-ssDNA filament assembly, and by helping UvsX to displace bound gp32 from ssDNA.

Gene ID:

uvy

Uniprot

P04537

Synonyms:

Recombination protein uvy, uvy

Immunogen:

Recombinant Enterobacteria phage T4 Recombination protein uvy protein (1-137AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A