

PACO34990

---

## Product Information

**Size:**

50ug

**Reactivity:**

Aspergillus niger

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Hydrolyzes a variety of simple alpha-D-galactoside as well as more complex molecules such as oligosaccharides and polysaccharides.

**Gene ID:**

aglB

**Uniprot**

Q9Y865

**Synonyms:**

Probable alpha-galactosidase B (EC 3.2.1.22) (Melibiase B), aglB

**Immunogen:**

Recombinant Aspergillus niger Probable alpha-galactosidase B protein (21-443AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A