

PACO34994

Product Information

Size:

50ug

Reactivity:

Aspergillus niger

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Hydrolyzes a variety of simple alpha-D-galactoside as well as more complex molecules such as oligosaccharides and polysaccharides.

Gene ID:

aglB

Uniprot

A2QEJ9

Synonyms:

Probable alpha-galactosidase B (EC 3.2.1.22) (Melibiase B), aglB

Immunogen:

Recombinant Aspergillus niger Probable alpha-galactosidase B protein (21-443AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A