

PACO35010

---

## Product Information

**Size:**

50ug

**Reactivity:**

Bradyrhizobium sp.

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Is an aliphatic amidase with a restricted substrate specificity, as it only hydrolyzes formamide.

**Gene ID:**

amiF

**Uniprot**

A4Z3G9

**Synonyms:**

Formamidase (EC 3.5.1.49) (Formamide amidohydrolase), amiF

**Immunogen:**

Recombinant Bradyrhizobium sp. Formamidase protein (1-337AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A