

mip Antibody, HRP conjugated



PACO35015

Product Information

Size:

50ug

Reactivity:

Legionella pneumophila

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Essential virulence factor associated with macrophage infectivity. Exhibits PPIase activity.

Gene ID:

mip

Uniprot

A5IGB8

Synonyms:

Outer membrane protein MIP (EC 5.2.1.8) (Macrophage infectivity potentiator) (Peptidyl-prolyl cis-trans isomerase) (PPIase) (Rotamase), mip

Immunogen:

Recombinant Legionella pneumophila Outer membrane protein MIP protein (21-233AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A