

PACO35018

Product Information

Size:

50ug

Reactivity:

Bacteroides vulgatus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Alpha-galactosidase. Removes both branched alpha-1,3-linked galactose residues of blood group B antigens and linear alpha-1,3-linked galactose structures.

Gene ID:

glaB1

Uniprot

A6KWM0

Synonyms:

Alpha-1,3-galactosidase B (EC 3.2.1. n1) (EC 3.2.1. n2) (Exo-alpha-galactosidase B1) (EC 3.2.1.22), glaB1

Immunogen:

Recombinant Bacteroides vulgatus Alpha-1,3-galactosidase B protein (19-592AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A