

COL6A6 Antibody



PACO35026

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Collagen VI acts as a cell-binding protein.

Gene ID:

COL6A6

Uniprot

A6NMZ7

Synonyms:

Collagen alpha-6(VI) chain, COL6A6


Immunogen:

Recombinant Human Collagen alpha-6(VI) chain protein (430-982AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO35026 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO35026 at dilution of 1:100.