

PACO35030

---

## Product Information

**Size:**

50ug

**Reactivity:**

Escherichia coli O9: H4

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Participates in the translocation of lipoproteins from the inner membrane to the outer membrane. Only forms a complex with a lipoprotein if the residue after the N-terminal Cys is not an aspartate (The Asp acts as a targeting signal to indicate that the lipoprotein should stay in the inner membrane).

**Gene ID:**

lolA

**Uniprot**

A7ZYJ5

**Synonyms:**

Outer-membrane lipoprotein carrier protein, lolA

**Immunogen:**

Recombinant Escherichia coli O9: H4 Outer-membrane lipoprotein carrier protein (22-203AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A