

# dsdA Antibody



PACO35042

---

## Product Information

**Size:**

50ug

**Reactivity:**

Salmonella heidelberg

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

dsdA

**Uniprot**

B4TA53

**Synonyms:**

D-serine dehydratase (EC 4.3.1.18) (D-serine deaminase) (DSD), dsdA

**Immunogen:**

Recombinant Salmonella heidelberg D-serine dehydratase protein (1-440AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A