

dsdA Antibody



PACO35046

Product Information

Size:

50ug

Reactivity:

Salmonella heidelberg

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

dsdA

Uniprot

B4TA53

Synonyms:

D-serine dehydratase (EC 4.3.1.18) (D-serine deaminase) (DSD), dsdA

Immunogen:

Recombinant Salmonella heidelberg D-serine dehydratase protein (1-440AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A