

PACO35050

Product Information

Size:

50ug

Reactivity:

Escherichia coli O127: H6

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Involved in protein export. Acts as a chaperone by maintaining the newly synthesized protein in an open conformation. Functions as a peptidyl-prolyl cis-trans isomerase.

Gene ID:

tig

Uniprot

B7UJQ9

Synonyms:

Trigger factor (TF) (EC 5.2.1.8) (PPlase), tig

Immunogen:

Recombinant Escherichia coli O127: H6 Trigger factor protein (1-432AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A