

PACO35062

Product Information

Size:

50ug

Reactivity:

Sulfolobus islandicus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Catalyzes the first reaction in the catabolism of the essential branched chain amino acid, leucine, isoleucine, and valine. May also function as a transporter of branched chain alpha-keto acid, .

Gene ID:

creN7

Uniprot

C3MPN0

Synonyms:

Chromatin protein Cren7, creN7

Immunogen:

Recombinant Sulfolobus islandicus Chromatin protein Cren7 protein (1-60AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A