

mgsA Antibody



PACO35066

Product Information

Size:

50ug

Reactivity:

Bacillus cereus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

mgsA

Uniprot

C1EN30

Synonyms:

Methylglyoxal synthase (MGS) (EC 4.2.3.3), mgsA

Immunogen:

Recombinant Bacillus cereus Methylglyoxal synthase protein (1-131AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A