

# blaNDM-1 Antibody, FITC conjugated



PACO35072

---

## Product Information

**Size:**

50ug

**Reactivity:**

Klebsiella pneumoniae

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Confers resistance to many beta-lactam antibiotics, including some carbapenems. Does not confer resistance to the polymixin colistin or the fluoroquinolone ciprofloxacin.

**Gene ID:**

blaNDM-1

**Uniprot**

C7C422

**Synonyms:**

Metallo-beta-lactamase type 2 (EC 3.5.2.6) (B2 metallo-beta-lactamase) (Beta-lactamase type II) (Metallo-beta-lactamase NDM-1) (Metallo-beta-lactamase type II) (New Delhi metallo-beta-lactamase-1) (NDM-1), blaNDM-1

**Immunogen:**

Recombinant Klebsiella pneumoniae Metallo-beta-lactamase type 2 protein (29-270AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A