

PACO35074

---

## Product Information

**Size:**

50ug

**Reactivity:**

Microbacterium testaceum

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Hydrolyzes short- and long-chain N-acyl homoserine lactones with or without 3-oxo substitution at C3. Has slight activity towards L-homoserine lactone, and no activity towards gamma-butyrolactone.

**Gene ID:**

aiiM

**Uniprot**

C6L862

**Synonyms:**

N-acyl homoserine lactonase (AHL-lactonase) (EC 3.1.1.81), aiiM

**Immunogen:**

Recombinant Microbacterium testaceum N-acyl homoserine lactonase protein (1-251AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A