

# gL Antibody, FITC conjugated



PACO35084

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human cytomegalovirus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

The heterodimer glycoprotein H-glycoprotein L is required for the fusion of viral and plasma membranes leading to virus entry into the host cell. Membrane fusion is mediated by the fusion machinery composed at least of gB and the heterodimer gH/gL.

**Gene ID:**

gL

**Uniprot**

F5HCH8

**Synonyms:**

Envelope glycoprotein L (gL), gL

**Immunogen:**

Recombinant Human cytomegalovirus Envelope glycoprotein L protein (31-278AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A