

RPN13 Antibody



PACO35086

Product Information

Size:

50ug

Reactivity:

Saccharomyces cerevisiae

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Component of the 19S cap proteasome complex which acts as a regulatory subunit of the 26S proteasome, involved in the ATP-dependent degradation of ubiquitinated proteins.

Gene ID:

RPN13

Uniprot

O13563

Synonyms:

26S proteasome regulatory subunit RPN13 (Proteasome non-ATPase subunit 13), RPN13, DAQ1

Immunogen:

Recombinant Saccharomyces cerevisiae 26S proteasome regulatory subunit RPN13 protein (2-156AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A