

PACO35122

---

## Product Information

**Size:**

50ug

**Reactivity:**

Acetobacter pasteurianus

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

This methylase recognizes the double-stranded sequence GTGCAC, causes specific methylation on C on both strands, and protects the DNA from cleavage by the ApaLI endonuclease.

**Gene ID:**

apaLIM

**Uniprot**

O52702

**Synonyms:**

Modification methylase ApaLI (M. ApaLI) (EC 2.1.1.37) (Cytosine-specific methyltransferase ApaLI), apaLIM

**Immunogen:**

Recombinant Acetobacter pasteurianus Modification methylase ApaLI protein (1-429AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A