

# B4GALT5 Antibody, FITC conjugated



PACO35156

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

Responsible for the synthesis of complex-type N-linked oligosaccharides in many glycoproteins as well as the carbohydrate moieties of glycolipids.

**Gene ID:**

B4GALT5

**Uniprot**

O43286

**Synonyms:**

Beta-1,4-galactosyltransferase 5 (Beta-1,4-GalTase 5) (Beta4Gal-T5) (b4Gal-T5) (EC 2.4.1.-) (Beta-1,4-GalT II) (UDP-Gal: beta-GlcNAc beta-1,4-galactosyltransferase 5) (UDP-galactose: beta-N-acetylglucosamine beta-1,4-galactosyltransferase 5), B4GALT5

**Immunogen:**

Recombinant Human Beta-1,4-galactosyltransferase 5 protein (101-388AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A