

B4GALT5 Antibody, Biotin conjugated



PACO35157

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Responsible for the synthesis of complex-type N-linked oligosaccharides in many glycoproteins as well as the carbohydrate moieties of glycolipids.

Gene ID:

B4GALT5

Uniprot

O43286

Synonyms:

Beta-1,4-galactosyltransferase 5 (Beta-1,4-GalTase 5) (Beta4Gal-T5) (b4Gal-T5) (EC 2.4.1.-) (Beta-1,4-GalT II) (UDP-Gal: beta-GlcNAc beta-1,4-galactosyltransferase 5) (UDP-galactose: beta-N-acetylglucosamine beta-1,4-galactosyltransferase 5), B4GALT5

Immunogen:

Recombinant Human Beta-1,4-galactosyltransferase 5 protein (101-388AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A