

NP Antibody



PACO35178

Product Information

Size:

50ug

Reactivity:

Influenza A virus

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

NP

Uniprot

O91743

Synonyms:

Nucleoprotein (Nucleocapsid protein) (Protein N), NP

Immunogen:

Recombinant Influenza A virus Nucleoprotein protein (1-498AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A