

SCRG1 Antibody, HRP conjugated



PACO35187

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Gene ID:

SCRG1

Uniprot

O75711

Synonyms:

Scrapie-responsive protein 1 (Scrapie-responsive gene 1 protein) (ScRG-1), SCRG1

Immunogen:

Recombinant Human Scrapie-responsive protein 1 protein (21-98AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A