

# SCRG1 Antibody, FITC conjugated



PACO35188

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

SCRG1

**Uniprot**

O75711

**Synonyms:**

Scrapie-responsive protein 1 (Scrapie-responsive gene 1 protein) (ScRG-1), SCRG1

**Immunogen:**

Recombinant Human Scrapie-responsive protein 1 protein (21-98AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A