

ABHD18 Antibody



PACO35242

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

ABHD18

Uniprot

Q0P651

Synonyms:

Protein ABHD18 (Alpha/beta hydrolase domain-containing protein 18) (Abhydrolase domain-containing protein 18), ABHD18, C4orf29


Immunogen:

Recombinant Human Protein ABHD18 protein (25-414AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO35242 at dilution of 1:100.