

RASGEF1B Antibody



PACO35246

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Guanine nucleotide exchange factor (GEF) with specificity for RAP2A, it doesn't seem to activate other Ras family proteins (in vitro).

Gene ID:

RASGEF1B

Uniprot

Q0VAM2

Synonyms:

Ras-GEF domain-containing family member 1B (GPI gamma-4), RASGEF1B, GPIG4


Immunogen:

Recombinant Human Ras-GEF domain-containing family member 1B protein (1-300AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO35246 at dilution of 1:100.