

ARL16 Antibody



PACO35258

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

ARL16

Uniprot

Q0P5N6

Synonyms:

ADP-ribosylation factor-like protein 16, ARL16


Immunogen:

Recombinant Human ADP-ribosylation factor-like protein 16 protein (1-197AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO35258 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO35258 at dilution of 1:100.