

# LY6D Antibody



PACO35298

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May act as a specification marker at earliest stage specification of lymphocytes between B- and T-cell development. Marks the earliest stage of B-cell specification.

**Gene ID:**

LY6D

**Uniprot**

Q14210

**Synonyms:**

Lymphocyte antigen 6D (Ly-6D) (E48 antigen), LY6D, E48

**Immunogen:**

Recombinant Human Lymphocyte antigen 6D protein (21-98AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A